



Animal Health Catalogue

Surgical clippers





 CareFusion

WET/DRY



Surgical clippers

- A single handle with easy to change clipping heads :
General use blade, Neuro blade for thickest hair,
SensiClip™ blade for sensitive or loose skin
- The cutting ability, closeness of cut and specific head types for each hair type allow to remove hair in a single pass
- 40° angled ergonomic handle creates a more intuitive use for seamless hair removal in a single pass
- Clipping closer to the skin without disrupting it, 0.21mm cut height
- Wet/dry capability
- The full submersibility of the clipper ensures thorough cleaning and disinfection.
- Made of stainless steel/ABS* polymer for resistance and toughness
- Lithium-ion battery technology facilitates longer run time and improves energy efficiency
- Simple assembly of lids and containers by special latching
- Can be used in wet or dry clipping conditions and is fully submersible* for enhanced disinfection, easy management and better efficiency in your practice
- Single use disposable blades



Description	Units/ Box	Reference
Rechargeable surgical clipper	1	5513E
Charging adapter - EU	1	5514E
Charging adapter - UK	1	5514U
Neuro blade	20	4412A
SensiClip® blade	20	4403A



Neuro blade



SensiClip blade



Learn more about our Surgical Clippers by scanning the QR Code with your smart phone.

BD Switzerland Sàrl, Terre Bonne Park – A4, Route de Crassier 17, 1262 Eysins, Switzerland
Tel: +41 21 556 30 00

[bd.com](https://www.bd.com)

© 2019 BD. BD and the BD Logo are trademarks of Becton, Dickinson and Company. CF04131

